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

サンプル	Sequence Entry	Accession No	%
<b>トウモロコシ</b> (学名: <i>Zea mays</i> ) 	<i>Zea mays</i> clone 1503689 mRNA sequence	EU955045	96.65
	<i>Zea mays</i> cultivar line X17 16S ribosomal RNA,16S-23S ribosomal RNA intergenic spacer, partial sequence	DQ369917	96.3
	<i>Zea mays</i> cultivar line T62 16S ribosomal RNA,16S-23S ribosomal RNA intergenic spacer, partial sequence	DQ369910	96.3
	<i>Zea mays</i> clone Contig488 mRNA sequence	BT016655	95.91
	<i>Zea mays</i> subsp. <i>mays</i> x <i>Zea perennis</i> 16S ribosomal RNA, 16S-23S ribosomal RNA intergenic spacer,partial sequence	DQ369908	95.17
<b>大豆</b> (学名: <i>Glycine max</i> ) 	<i>Glycine max</i> 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence; and 26S rribosomal RNA gene, partial sequence	FJ609734	100
	<i>Glycine max</i> cultivar Hwangkeum 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, complete sequence; and 5.8S ribosomal RNA gene,partial sequence	EU118312	100
	<i>Glycine soja</i> cultivar IT182932 18S ribosomal RNA ene,partial sequence; internal transcribed spacer 1, complete sequence; and 5.8S ribosomal RNA gene, partial sequence	EU118313	100
	<i>Glycine max</i> voucher PS1684MT01 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2,complete sequence; and 28S ribosomal RNA gene, partial sequence	FJ980442	100
	<i>Glycine max</i> isolate 48Hidatsa internal transcribed spacer 1, partial sequence; 5.8S ribosomal RNA gene, complete sequence; and internal transcribed space 2, partial sequence	AF144654	99.47
<b>キウイ</b> (学名: <i>Actinidia deliciosa</i> あるいは <i>Actinidia chinensis</i> ) 	<i>Actinidia deliciosa</i> genes for 18S ribosomal RNA, internal transcribed spacer 1, 5.8S ribosomal RNA, partial and complete sequence	AB253775	99.66
	<i>Actinidia henanensis</i> internal transcribed spacer 1, partial sequence; 5.8S ribosomal RNA gene,complete sequence; and internal transcribed spacer 2, partial sequence	AF323841	99.62
	<i>Actinidia deliciosa</i> clone 1 internal transcribed spacer 1, partial sequence; 5.8S ribosomal RNA gene, complete sequence; and internal transcribed spacer 2, partial sequence	AF323830	99.62
	<i>Actinidia callosa</i> var. <i>strigillosa</i> clone 20 internal transcribed spacer 1, partial sequence; 5.8S ribosomal RNA gene, complete sequence; and internal transcribed spacer 2, partial sequence	AF323829	99.62
	<i>Actinidia callosa</i> var. <i>strigillosa</i> clone 10 internal transcribed spacer 1, partial sequence; 5.8S ribosomal RNA gene, complete sequence; and internal transcribed spacer 2, partial sequence	AF323827	99.61

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サンプル	Sequence Entry	Accession No	%
トマト ( <i>Solanum lycopersicum</i> ) 	Solanum lycopersicum genes for 18S rRNA, ITS1, 5.8S rRNA, ITS2, 28S rRNA, partial and complete sequence, cultivar: Momotaro	AB373816.1	100%
	Solanum lycopersicum genes for 18S rRNA, ITS1, 5.8S rRNA, ITS2, 28S rRNA, partial and complete sequence, note: jtomate creollo	AB373815.1	100%
	Solanum lycopersicum var. cerasiforme genes for 18S rRNA, ITS1, 5.8S rRNA, ITS2, 28S rRNA, partial and complete sequence	AB373814.1	100%
	Solanum pimpinellifolium genes for 18S rRNA, ITS1, 5.8S rRNA, ITS2, 28S rRNA, partial and complete sequence	AB373813.1	100%
	Lycopersicon esculentum internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence	AF244747.1	100%
リンゴ (学名: <i>Malus pumila</i> ) 	Malus domestica cultivar Ashmeads Kernal internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence	AF186480.1	97%
	Malus domestica cultivar Reinette Simerenko internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence	AF186478.1	97%
	Malus sieversii isolate sieversiiJ5 internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence	AF186490.1	97%
	Malus sieversii isolate sieversiiJ8 internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence	AF186493.1	97%
	Malus sieversii isolate sieversiiJ7 internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence	AF186492.1	97%
グレープフルーツ (学名: <i>Citrus x paradisi</i> ) 	Citrus x paradisi bio-material SOFRI:STGBC603 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence; and 28S ribosomal RNA gene, partial sequence	FJ641956.1	99%
	Citrus maxima bio-material BIRDI:SH-TKD5 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence; and 28S ribosomal RNA gene, partial sequence	FJ641954.1	99%
	Citrus maxima bio-material BIRDI:SH-P36 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence; and 28S ribosomal RNA gene, partial sequence	FJ641953.1	99%
	Citrus maxima bio-material BIRDI:SH-P9 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence; and 28S ribosomal RNA gene, partial sequence	FJ641952.1	99%
	Citrus maxima bio-material FRDC:SH-P47 18S ribosomal RNA gene, partial sequence; internal transcribed spacer 1, 5.8S ribosomal RNA gene, and internal transcribed spacer 2, complete sequence; and 28S ribosomal RNA gene, partial sequence	FJ641951.1	99%

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